



Powder Free Vinyl Gloves 100

	Small	Medium	Large	X Large
Product Code	67045	70642	31376	35664



Material Specification	Clear Vinyl (Polyvinyl Chloride - PVC)	
Manufacturing Process Type	Straight Vinyl Dip With Polyurethene Coating	
Colour	Clear	
Texture & Grip	Non Textured Surface	
Former Type	Non Textured Former	

BS EN 455 - 2 Phy	sical properties
--------------------------	------------------

	Small	Medium	Large	X Large
Palm Width +/- 5mm	85mm	95mm	105mm	115mm
Single Wall Thickness				
Finger	0.08mm	0.08mm	0.08mm	0.08mm
Palm	0.05mm	0.05mm	0.05mm	0.05mm
Cuff	0.04mm	0.04mm	0.04mm	0.04mm
Weight Of Glove	4g	5g	6g	7g
Tensile Strength +/- 5mpa	17mpa	17mpa	17mpa	17mpa

BS EN 455 - 3	Biological Evaluation
---------------	-----------------------

AQL 1.5 - Visual Defects Aceptance Criteria				
Major Visual Defects	Insect	Torn	Touching	
	Former knocking	Visible hole		
Minor Visual Defects	Thin spots	Stain / Dirt ≤ 2.5mm ²	Sticky	
AQL 2.5	Mix size	Mix type	Donning	
	Mixed colour	Discolouration	Dirt / Stain ≥ 2.5mm ²	
	Sticky	Wet Look	Lumps ≥ 2.5mm ²	
	Hair / Foreign matter			
AQL 4.0	Wet look	Lump ≤ 2.5mm ²	Poor beading	
	Double dipped			

Quality Standards	Manufacturing QMS	BS EN ISO 9001 / BS EN ISO 13485
	Glove	BS EN 455 Parts 1, 2 & 3
	Sampling	ISO 2859-1
	Medical Device Directive Requirement	93/42/EEC Class 1
	Shelf Life	3 years