

# AURELIA® Ignite®

Examination Gloves



97887	(M)	ø 100	10 x 100
97888	(L)	ø 100	10 x 100
97889	(XL)	ø 90	10 x 90
97880	(XXL)	ø 90	10 x 90



Nitrile



Powder Free



Orange



8.4 +/- 0.2 g



Diamond-Textured




5 Years

Premium quality heavy-duty Orange Nitrile examination gloves with an excellent diamond grip finish which offers maximum grip, increased tensile strength, optimum protection in tough environments and superior tactile sensitivity.

- 7mil Powder-Free Nitrile
- Non-Sterile
- Beaded cuff
- Ambidextrous
- Latex Free
- AQL 1.5

## Sizes & Dimensions



M	95 ± 10 mm	min. 240 mm
L	110 ± 10 mm	min. 240 mm
XL	120 ± 10 mm	min. 240 mm
XXL	≥ 120 mm	min. 240 mm

## Single wall thickness



- ① min. 0.25 mm
- ② min. 0.20 mm
- ③ min. 0.09 mm

## Physical properties

	Unaged	Aged
Tensile Strength	min. 14.0 MPa	min. 14.0 MPa
Elongation	min. 500 %	min. 400 %
Strength at Break	min. 6N	min. 6N

## Quality Signs & Certificates

89/686/EEC, 93/42/EWG

ISO 9001:2008, ISO 13485, CMDCAS 13485 (SGS)

EN 455 1-4

EN 420 : 2003 + A 1 : 2009

EN 374 2-3, EN 388

EN 1186: Food Contact Approved (SGS)

