

Approved Manufacturers

Clinell Universal Wipes have been approved for use by the following manufacturers:

BRAND	EQUIPMENT
1st Call Mobility	Mattresses
Abbott	FreeGo Pump
Abecca	Positioning Supports & Mattresses
Active Key GmbH & Co. KG	All Active Key MedicalKey™ products, including: AK-C4112, AK-C4412, AK-C7012 & AK-C7412
AGFA	DX-D Products
Agile Medical	All Products
AKV International	Upholstery Vinyl
Albyn Medical	Hermes Urodynamic (do not use wipes on monitor)
Align Technology	iTero Element
All Modul	Trolley Range
Arcomed	Pumps
ArjoHuntleigh	Mattresses Reliant IS2 Cover (this is a mattress cover)
ASep Healthcare	Tournistretch
BD (Bard)	Site-Rite® 8 Ultrasound System and Probe CareFusion Jaeger Spirometry Cabinet Alaris Syringe Pump Alaris GW Volumetric Pump
Bedfont	Smokerlyzer, ToxCO (not suitable with Gastrolyzer)
Brandon Medical	Astralite Treatment Lights Astramax Minor Surgical Lights Quasar Elite Operating Theatre Lights

Approved Manufacturers

Clinell Universal Wipes have been approved for use by the following manufacturers:

BRAND	EQUIPMENT
Braun	Thermoscans
Carrflex	Mattress Covers, Upholstery Tops & Associated Equipment
Covidien	Genius 2 Tympanic Thermometer
Dartex Coating	Mattresses (with polyurethane face)
De Smit Medical	CubeScan Bladder Scanners
Direct Healthcare	Dartex Mattress Covers
Echosens	FibroScan
Elekta	IBeam Evo Couch Top
Fresenius Kabi - Agilia range	VP MC Volumetric Infusion Pump SP MC Syringe Infusion Pump SP TIVA TCI Syringe Infusion Pump SP PCA Syringe Infusion Pump
Fritz Stephen	Pediatric Ventilation Systems
Fujifilm Sonosite	Ultrasound Equipment (Edge, NanoMaxx, MicroMaxx, TITAN) iViz Devices and C60v, L25v, L38v & P21v Probe

Approved Manufacturers

Clinell Universal Wipes have been approved for use by the following manufacturers:

BRAND	EQUIPMENT
GE Healthcare	<p>Ultrasound Transducers (please email info@gamahealthcare.com for full list of compatible transducers) Probe Holders, User interface, Touch Panel, OLED Monitor Display and the Housing of Voluson Expert BT15, BT16, BT17, BT18 and BT19 Consoles. Note: LCD Monitor Display can only be cleaned with IPA (70%) solution Voluson Ultrasound Consoles Expert Series (E6, E8 and E10) Signature Series (S8 and S10) Performance Series (P8) MAC VU360™ Resting ECG</p>
Haines Medical Australia	<p>Haines MediCurtains (full range) Haines Wipeclean Waterproof Pillows Haines Reusable Tourniquets Haines PVC Wipedown Mattress Covers Haines PVC Mackintosh Rolls Haines Reusable Pillow Cases</p>
Hospedia	<p>Patient Entertainment Unit</p>
Hospira	<p>Sapphire Pump</p>
Howard Wright	<p>Mattresses</p>

Approved Manufacturers

Clinell Universal Wipes have been approved for use by the following manufacturers:

BRAND	EQUIPMENT
HP	IPAQ HP EliteBook 840 G6 Healthcare Edition HP EliteOne 800 G5 Healthcare Edition All-in-One HP Healthcare Edition HC270cr Clinical Review Display HP Healthcare Edition HC241 Clinical Review Monitor HP Healthcare Edition HC241p Clinical Review Monitor HP Healthcare Edition HC271 Clinical Review Monitor HP Healthcare Edition HC271p Clinical Review Monitor
Huntleigh Healthcare	Foetal Monitors
International Light Technologies	IL1350
Karomed	Mattresses
Konica Minolta	SONOVISTA FX Ultrasound Device (premium edition)
Laerdal	Suction Units
MePACS	Duress Alarms
Mindray	N Series Monitor

Approved Manufacturers

Clinell Universal Wipes have been approved for use by the following manufacturers:

BRAND	EQUIPMENT
Modsel	All Product Ranges
OES Medical	Pressure Monitor and UAM Anaesthetic Machine
Orfit Industries	All Related Types of AIO 3.0 Medium Density Blue Head Support
OxyLitre	All Equipment
Panaz	Cadet Furniture
Parity Medical	Mobile Computer Carts
Park House	Healthcare Mattress Covers
Phantom Laboratory	Catphan 500
Philips	Intellivue Monitor (remaining residue must be wiped off) Ultrasound Transducers (please email info@gamahealthcare.com for full list of compatible transducers)
Plinth Medical	Vinyl Couches
Recticel	Mattresses
Seca	All Equipment
Sepura	Sepura TETRA Products
Sidhil	IQ Bed Frame
Siemens	Mammomat Inspiration
Smiths Medical	Pressure Infuser
Toshiba	Ultrasound Equipment
Tough-PAC	iPad Covers
Triton Electronics Systems Ltd	Multigas Analyser-AMG-06

Approved Manufacturers

Clinell Universal Wipes have been approved for use by the following manufacturers:

BRAND	EQUIPMENT
Universal Ferrule Remover	Ferrule Remover
Verathon	GlideScope Cables BladderScan Prime
Visionflex	Thinklabs One Digital Stethoscope
Welch Allyn	Connex Spot Monitor Caretemp Thermometers Pro6000 Thermometers Pocket Plus LED Diagnostic Set - Ophthalmoscope
Zebra	Zebra Healthcare Mobile Computers
Zoll	AED