



AURELIA® BOLD MAX

6mil, Nitrile Powder Free Examination Gloves



AURELIA® GLOVES

A SUPERMAX HEALTHCARE BRAND


sales@supermax.co.uk

+441733 615215

Supermax Healthcare Ltd,
Unit 12-16 Titan Drive
Fengate, Peterborough
PE1 5XN

aureliagloves.co.uk

Regulations	PPE Reugulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms. BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. EN 1184: 1997. Materials and articles in contact with foodstuffs.



Sizes & Dimensions	Single wall thickness	Physical Properties												
M 80+10 mm L 95+10 mm XL 110+10 mm XXL > 110 mm	min. 240 mm min. 240 mm min. 240 mm min. 240 mm	 ① min. 0.20 mm ② min. 0.16 mm ③ min. 0.08 mm												
		<table border="1"><tr><th></th><th>Unaged</th><th>Aged</th></tr><tr><td>Tensile Strength</td><td>Min. 14.0 MPa</td><td>Min. 14.0 MPa</td></tr><tr><td>Elongation</td><td>Min. 500%</td><td>Min. 400%</td></tr><tr><td>Strength at Break</td><td>Min. 6N</td><td>Min. 6N</td></tr></table>		Unaged	Aged	Tensile Strength	Min. 14.0 MPa	Min. 14.0 MPa	Elongation	Min. 500%	Min. 400%	Strength at Break	Min. 6N	Min. 6N
	Unaged	Aged												
Tensile Strength	Min. 14.0 MPa	Min. 14.0 MPa												
Elongation	Min. 500%	Min. 400%												
Strength at Break	Min. 6N	Min. 6N												



Nitrile



Powder Free



5 Years



Diamond



7.5g +/- 0.2g



Black

