

Product specification sheet

Revision dated 20 September 2019

1.	Product code	GS008*-A									
2.	Product description	Unicare Blue Vinyl									
3.	Material	Vinyl									
4.	Inspection parameters	Criteria		Inspection level	AQL	Standard					
Dimensions		S-2	4.0	EN455-2							
Physical properties		S-2	4.0	EN455-2							
Water tight test (1,000ml)		G-1	1.5	EN455-1							
Visual defects		G-1	2.5	In-house							
5.	Glove dimensions	Dimensions	Size (mm)								
XS			S	M	L	XL					
Product code		-	GS0082-A	GS0083-A	GS0084-A	GS0085-A					
Palm width		-	80 ± 10	95 ± 10	110 ± 10	≥ 110					
Length		-	≥ 240	≥ 240	≥ 240	≥ 240					
6.	Glove thickness	Point of thickness measurement			Nominal single wall (mm)						
Cuff (at beaded end)			0.05								
Palm (at centre of palm)			0.06								
Fingertip (at fingertip)			0.07								
7.	Physical properties	Parameters		Before aging	After aging						
Median force at break (N)		Min 3.6		Min 3.6							
Tensile strength (MPa)		Min 9		Min 9							
Ultimate elongation (%)		Min 300		Min 300							
8.	Proteins	None									
9.	Powder residue	Maximum 2.0 mg / glove									
10.	Cuff type	Beaded									
11.	Design	Ambidextrous									
12.	Surface	Smooth									
13.	Colour	Blue									
14.	Quality assurance	Manufactured under ISO 9001 and ISO 13485									
15.	Weight (size medium)	4.5g									
16.	Standards compliant with	EN455-1	EN455-2	EN455-3	EN455-4	EN420	EN374-1	EN374-2	EN374-4	EN374-5	Food approved
		✓	✓	✓	✓	✓					✓

The above data represents the product norm, which may be changed without notice and is intended as a guide only.