

Product specification sheet

Revision dated 20 September 2019

1.	Product code	GS012*									
2.	Product description	Unicare Green Vinyl									
3.	Material	Vinyl									
4.	Inspection parameters	Criteria		Inspection level		AQL		Standard			
		Dimensions		S-2		4.0		EN455-2			
		Physical properties		S-2		4.0		EN455-2			
		Water tight test (1,000ml)		G-1		4.0		EN455-1			
		Visual defects		G-1		2.5		In-house			
5.	Glove dimensions	Dimensions		Size (mm)							
				XS	S	M	L	XL			
		Product code		-	GS0122	GS0123	GS0124	GS0125			
		Palm width		-	80 ± 10	95 ± 10	110 ± 10	≥ 110			
		Length		-	≥ 240	≥ 240	≥ 240	≥ 240			
6.	Glove thickness	Point of thickness measurement					Nominal single wall (mm)				
		Cuff (at beaded end)					0.05				
		Palm (at centre of palm)					0.06				
		Fingertip (at fingertip)					0.07				
7.	Physical properties	Parameters			Before aging			After aging			
		Median force at break (N)			N/A			N/A			
		Tensile strength (MPa)			Min 9			Min 9			
		Ultimate elongation (%)			N/A			N/A			
8.	Proteins	None									
9.	Powder residue	Maximum 2.0 mg / glove									
10.	Cuff type	Beaded									
11.	Design	Ambidextrous									
12.	Surface	Smooth									
13.	Colour	Green									
14.	Quality assurance	Manufactured under ISO 9001 and ISO 13485									
15.	Weight (size medium)	4.5g									
16.	Standards compliant with	EN455-1	EN455-2	EN455-3	EN455-4	EN420	EN374-1	EN374-2	EN374-4	EN374-5	Food approved
						✓					✓

The above data represents the product norm, which may be changed without notice and is intended as a guide only.