

# PRODUCT SPECIFICATION



## ØRN Harpy Vinyl PF Blue Examination Gloves

**PRODUCT CODE:**

OH1521

**COLOUR:**

Blue

**SIZES:**

S, M, L, XL

**DESCRIPTION:**

- Non-Sterile
- Powder Free
- Ambidextrous
- Compliant with Food Handling Requirements
- Latex-Free
- Dop or Dehp Free
- Higher Tensile Strength and More Tactile Sensitivity
- Beaded Cuff for Added Strength and Easier Wearing.

**PACKAGING:**

100 Gloves Per Box  
10 Boxes Per Case  
77 Cartons Per Pallet



**CERTIFICATIONS:**

- CE: CAT III PPE
- UKCA
- EN 455-1:2020
- EN 455-2:2015
- EN 455-3:2015
- EN 455-4:2009
- EN ISO 374-1:2020
- EN ISO 374-2:2019
- EN ISO 374-4:2019
- EN ISO 374-5: 2016
- EN ISO 16523-1:2015
- EN 21420
- AQL 1.5
- Foodsafe

**MANUFACTURING STANDARDS:**

- ISO 9001:2015
- ISO 13485:2016

**THICKNESS:**

3mil  
FINGER (mm): 0.08±0.03  
PALM (mm): 0.08±0.03  
CUFF (mm): 0.06±0.03

**LENGTH:**

240mm

**WEIGHT:**

4.5g +/- 0.3g

**SHELF LIFE:**

Up to 5 Years – EN455-4

**ORIGIN:**

China

**COMMODITY CODE:**

4015120000

