



### PRODUCT INFORMATION

Product Colour	Pink
Size Range	S-L
Palm Thickness	0.07 mm
Typical Weight	3.50 g
Shelf Life	5 Years

### FEATURES

**Chemical Resistant:** Achieves type B chemical resistance according to EN ISO 374-1, for protection against at least 3 different chemicals

**EN ISO 374-5 Protection:** Conforms with the requirements of EN ISO 374-5, for protection against bacteria, fungi and virus

**Comfort:** Nitrile technology provides an exceptionally soft and pliable material

**Grip:** Textured finish to the fingers to enhance grip and handling

**Food Safe:** Suitable for food contact; all food types. Please refer to our online Food Compatibility Chart for full details

**Sensitivity:** Latex free, reducing the risk of hypersensitivity and allergic reactions to natural rubber proteins

**Powder Content:** Powder Free, reducing the risk of dust contamination

**Silicone Free:** Silicone residues from gloves can interfere with painting, coating, or adhesive bonding processes, potentially causing defects. These silicone-free gloves help minimise that risk, ensuring clean application surfaces.

### SECTORS

Hair and Beauty • Janitorial • Medical, Vet, and Dental • Pharmaceutical and Laboratory • Retail and Household



### STANDARDS

Tested to European standards EN ISO 21420, EN ISO 374-1, EN ISO 374-4 & EN ISO 374-5

AQL 1.5

EU Certified

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425

Tested in compliance with Commission Regulation (EU) No 10/2011 for Food contact

**Compliant with European Standards EN455:** Medical Gloves for Single Use parts 1,2,3&4

Compliant with Medical Devices Regulation (EU) 2017/745 Class 1



MDR - Class I  
PPE CAT III

2777



CLASS 1  
MEDICAL  
DEVICE



EN ISO 374-1:  
2016+A1:  
2018 Type B

KPT



EN ISO 374-5:  
2016

VIRUS



For further technical information, please contact Customer Services on 01733 230 230 or the Technical team on 0333 320 8555



### CAUTION

Contains accelerator chemicals which may cause type IV skin irritation and contact dermatitis in sensitised individuals.

### PACKING

100 gloves per dispenser  
10 dispensers per case

### STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. You are advised to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear.

### DISPOSAL INSTRUCTIONS

Dispose of responsibly and in line with local guidance.

Size	Inner Packaging Barcode
Small	5024951031499
Medium	5024951031505
Large	5024951031512

Inner Packaging	
Inner dimensions (LxWxH)	220x120x61mm

Size	Outer Packaging Barcode
Small	5024951031468
Medium	5024951031475
Large	5024951031482

Outer Packaging	
Outer Dimensions (LxWxH)	315x230x248mm

### PLEASE NOTE:

Weight does not form part of the specification and is subject to change. Palm thicknesses is approximate and based on an agreed manufacturing tolerance. All thickness and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polycy Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 01733 230 230 or the Technical team on 0333 320 8555