



PRODUCT INFORMATION

Product Colour	Black
Size Range	S-XL
Fingertip Thickness	0.10 mm
Palm Thickness	0.08 mm
Typical Weight	4.00 g
Shelf Life	5 Years

FEATURES

Chemical Resistant: Achieves type C chemical resistance according to EN ISO 374-1, for chemical protection against at least 1 chemical

EN ISO 374-5 Protection: Conforms with the requirements of EN ISO 374-5, for protection against bacteria & fungi

Comfort: Nitrile technology provides an exceptionally soft and pliable material

Grip: Textured finish to the fingers to enhance grip and handling

Food Safe: Suitable for food contact; all food types. Please refer to our online Food Compatibility Chart for full details

Sensitivity: Latex free, reducing the risk of hypersensitivity and allergic reactions to natural rubber proteins

Powder Content: Powder Free, reducing the risk of dust contamination

Silicone Free: Silicone residues from gloves can interfere with painting, coating, or adhesive bonding processes, potentially causing defects. These silicone-free gloves help minimise that risk, ensuring clean application surfaces.

SECTORS

Automotive • Construction and Utility • Engineering and Manufacturing • Food Manufacturing and Preparation • Hair and Beauty • Janitorial • Medical, Vet, and Dental



STANDARDS

Tested to European standards EN ISO 21420, EN ISO 374-1, EN ISO 374-4 & EN ISO 374-5

AQL 1.5

EU Certified

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425

Compliant with Medical Devices Regulation (EU) 2017/745 Class 1

Tested in compliance with Commission Regulation (EU) No 10/2011 for Food contact

Compliant with European Standards EN455: Medical Gloves for Single Use parts 1,2,3&4




MDR - Class I
PPE CAT III



CLASS 1
MEDICAL
DEVICE



EN ISO 374-1:
2016+A1:
2018 Type C



EN ISO 374-5:
2016



2777

VIRUS



For further technical information, please contact Customer Services on 01733 230 230 or the Technical team on 0333 320 8555



CAUTION

Contains accelerator chemicals which may cause type IV skin irritation and contact dermatitis in sensitised individuals.

PACKING

100 gloves per dispenser
10 dispensers per case

STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. It is recommended to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear and tear.

DISPOSAL INSTRUCTIONS

Dispose of the product responsibly in accordance with local regulations.

Size	Inner Packaging Barcode
Small	5024951031567
Medium	5024951031574
Large	5024951031581
XL	5024951031598

Size	Outer Packaging Barcode
Small	5024951031529
Medium	5024951031536
Large	5024951031543
XL	5024951031550

For further technical information, please contact Customer Services on 01733 230 230 or the Technical team on 0333 320 8555