

PRODUCT SPECIFICATION



ØRN Golden Diamond Grip Nitrile Orange Examination Gloves

PRODUCT CODE:

OH1508

COLOUR:

Orange

SIZES:

M, L, XL, 2XL

DESCRIPTION:

- Non-Sterile
- Powder Free
- Ambidextrous
- Latex-Free
- Low In Residual Chemical
- Tear Resistant Providing Greater Durability Compared to Latex or Vinyl.
- Beaded Cuff for Added Strength and Easier Wearing.

PACKAGING:

50 Gloves Per Box
10 Boxes Per Carton
77 Cartons Per Pallet



CERTIFICATIONS:

- CE
- UKCA
- EN 455-1:2020
- EN 455-2:2015
- EN 455-3:2015
- EN ISO 374-1:2016+A1:2018/TYPE B/ KPTLO
- EN 16523-1:2015+A1:2018
- EN ISO 374-2:2019
- EN ISO 374-4:2019
- EN ISO 374-5:2016
- EN 21420:2020
- AQL 1.5

MANUFACTURING STANDARDS:

- ISO 9001:2015
- ISO 13485:2016

THICKNESS:

7.9mil
FINGER (mm): 0.20±0.03
PALM (mm): 0.20±0.03
CUFF (mm): 0.10±0.03

LENGTH:

240mm

WEIGHT:

8.0g +/- 0.3g

SHELF LIFE:

Up to 5 Years – EN455-4

ORIGIN:

China

COMMODITY CODE:

4015120000

DIAMOND
GRIP

