

PRODUCT SPECIFICATION



ØRN Cassin Clear Vinyl Examination Gloves

PRODUCT CODE:

OH1520

COLOUR:

Clear

SIZES:

S, M, L, XL, 2XL

DESCRIPTION:

- Non-Sterile
- Powder Free
- Beaded Cuff
- Ambidextrous
- Latex Free
- Dop or Dehp Free
- Higher Tensile Strength and More Tactile Sensitivity
- Compliant with Food Handling Requirements

PACKAGING:

100 Gloves Per Box
10 Boxes Per Carton
77 Cartons Per Pallet



CERTIFICATIONS:

- CE
- UKCA
- EN 455-1:2020
- EN 455-2:2015
- EN 455-3:2015
- EN ISO 374-1:2016: 2018/ TYPE B/ KPT
- EN 455-4:2009
- EN ISO 374-2:2019
- EN ISO 374-4:2019
- EN ISO 374-5: 2016
- EN ISO 16523-1:2015
- AQL 1.5
- Foodsafe

MANUFACTURING STANDARDS:

- ISO 9001:2015
- ISO 13485:2016

THICKNESS:

3mil
FINGER (mm): 0.08±0.03
PALM (mm): 0.08±0.03
CUFF (mm): 0.06±0.03

LENGTH:

240mm

WEIGHT:

4.5g +/- 0.3g

SHELF LIFE:

Up to 5 Years – EN455-4

ORIGIN:

China

COMMODITY CODE:

4015120000

