

Black Nitrile

Gloves Powder Free
YG39

Sizes

XS - XL

YG39 Datasheet Powder Free Black Nitrile Gloves



Glove Specifications - EN455

Size	Length mm	Palm Width mm	Thickness		
			Cuff	Palm	
XS	240	75 ± 10	Min 0.05	Min 0.05	Min 0.05
S	240	80 ± 10	Min 0.05	Min 0.05	Min 0.05
M	240	95 ± 10	Min 0.05	Min 0.05	Min 0.05
L	240	110 ± 10	Min 0.05	Min 0.05	Min 0.05
XL	240	≥110	Min 0.05	Min 0.05	Min 0.05

Performance Requirements

Test	Related Defects	Inspection Level	AQL
Water Leak at Final Release	Holes	G-1 Single Sampling Plan	1.5

Manufacturing Standards

Manufacturing Standards	In compliance with ISO 9001 and ISO 13485
Presentation	Standard: 100 gloves per dispenser box; 10 dispenser boxes per Carton
Glove Care	Should not store under direct sunlight or at high heat and humid conditions. Single use only.

General Details

Material	Acrylonitrile Butadiene (Nitrile) non sterile
Surface Finish	Powder Free
Colour	Black
Geometry	Ambidextrous
Texture	FingerTextured
Cuff and Finishing	Beaded



AQL 1.5

EN455 Parts 1, 2, 3 & 4

EN374

ASTM F1671-97

Class 1 Medical Device

MDD 93/42/EEC

